

# THERABATH PROFESSIONAL PARAFFIN WAX BATH - MODEL TB6

For over 90 years paraffin moist heat therapy has been used to relieve pain, inflammation and stiffness caused by arthritis. A form of thermotherapy, paraffin treatments effectively increase blood flow to aid in rehabilitation, or to warm an area of the body prior to other therapies or exercise.

The Therabath Professional Paraffin Bath is USA made to Class II medical device standards. Therabath units are hand-assembled ensuring that each and every one leaves the factory durable, safe, and reliable.

## MODEL TB6 FEATURES:

### DEEP IMMERSION 6 - 9 LB. WAX CAPACITY

up to 25% deeper than previous models

### ANODIZED ALUMINUM TANK

lightweight, durable, and energy efficient

### ALWAYS SAFE TEMPERATURE RANGE

even during the initial melt cycle!

### INTEGRATED HANDLES

ease transport in the home, clinic or office



Lid with Integrated Handle



5-Year Limited Warranty!

Includes 6 lbs (2.72 kg.) of Therabath Professional Refill Paraffin wax beads in your choice of 12 scents, or hypoallergenic ScentFree. Therabath Paraffin Wax is blended for optimal performance in the TB6.

Made in the USA and backed by a 5 Year Limited Warranty\*!

## PRODUCT SPECIFICATIONS

- 9 lb (4.08 kg.) Wax Capacity
- Aluminum Powder Coated Tank
- Unit Dimensions: 14.3" L x 8.2" W x 7.5" H
- Power: 100-125V, 1.04-1.55A, 50/60Hz
- Heater: 103-194 watts
- FDA 60601 Class II Medical Device Standards
- ISO 13485 Quality Standards



Integrated Handles for Portability & Rubber Feet for Added Stability



Removable IEC Electrical Cord

\*Contact us for full warranty details.

**therabath**®  
PROFESSIONAL THERMOTHERAPY



Contact us to Order!  
651.604.8480 | [therabath@wrmed.com](mailto:therabath@wrmed.com)

Copyright © 2020. Therabath is a registered trademark of WR Medical Electronics Co. All rights reserved. MR139, Rev. 01/25.