



ADJUSTABLE TEMPERATURE PARAFFIN WAX BATH - MODEL TB9

Item # 2275 TB9 - Suggested Price: \$383.50

The Therabath Model TB9 Adjustable Temperature Paraffin Bath features a **NEW SAFE QUICK MELT** setting, melting the paraffin in half the time of our Model TB6, but in an **always-safe temperature range!** No need to wait for the wax temperature to come back down before safely dipping!

NEW QUICK MELT MODE

Safely Melts 6lbs (2.72 kg.) of Paraffin in Just 3-4 hrs!



Includes 6 lbs (2.72 kg.) of Therabath Professional Refill Paraffin (ScentFree) wax blended for optimal performance in the TB9.

Made in the USA and backed by a Limited Lifetime Warranty*!

	Switch Position	Temperature Range	Melt Time*
	Quick Melt	131 - 142°F (55 - 61°C)	3-4 hrs.
	Normal	126 - 134°F (52 - 56°C)	6-8 hrs.

* Melt times are approximate; based on 6lbs (2.72 kg.) of paraffin wax.

PRODUCT SPECIFICATIONS

- 9 lb (4.08 kg.) Wax Capacity
- Anodized Aluminum Tank
- Unit Dimensions: 14.3" L x 8.2" W x 7.5" H
- Power: 100-125V, 1.04-1.55A, 50/60Hz
- Heater: 103-194 watts
- FDA 60601 Class II Medical Device Standards
- ISO 13485 Quality Standards

ETL CLASSIFIED



Safe for Operation in Both Quick Melt and Normal Modes!



Lid with Integrated Handle



Integrated Handles for Portability & Rubber Feet for Added Stability



Removable IEC Electrical Cord

*Contact us for full warranty details.



Contact us to Order!
651.604.8480 | therabath@wrmed.com

Copyright © 2020. Therabath is a registered trademark of WR Medical Electronics Co. All rights reserved. MR136, Rev. 06/20.